



Sample Preparation Products

Application Note: H-CX702

Extraction of Cocaine from Hair

Introduction

This Application Note describes a sample preparation method for the extraction of Cocaine from Hair. This method uses TELOS™ H-CX 130mg/3ml SPE Columns.

Sample Considerations

Analytes

Cocaine, Ecgonine methyl ester, Cocaethylene, Benzoyllecgonine.

Matrix

Hair.

Sample Preparation

Columns

TELOS H-CX 130mg/3ml SPE Columns, part number 500-130M-003T. For larger sample volumes, other sorbent mass/reservoir volume options are available.

Sample Pre-Treatment

Take 100mg of hair and treat with the following washing steps:

- 3 x 4ml of distilled water
- 3 x 4ml of acetone
- 3 x 4ml of ethanol
- 3 x 4ml of dichloromethane

N.B. If an initial positive identification test is required, it is possible to test the hair washings using a high sensitivity immunoassay. For this, combine the ethanol and dichloromethane washes, evaporate carefully to dryness (<40°C) under nitrogen and reconstitute as directed by the chosen test.


Dry the washed sample, cut into segments of 1mm and pulverise in a ball mill or similar. Place 50mg of pulverised sample in a 100 x 16mm tube, add 0.1M HCl (1ml) and heat at 37°C overnight. Add 0.1M phosphate buffer pH 6 (500>1) and homogenise in an agitator. If necessary, adjust the pH to 6.0 ± 0.5 with 100mM phosphate buffer pH 6.0.

Column Conditioning

Condition the column with methanol (3ml) at a flow rate of 4ml/min.

Column Equilibration

Equilibrate the column with water (3ml) at a flow rate of 4ml/min, followed by 0.1M phosphate buffer pH 6.0 (3ml) at the same flow rate.





Extraction of Cocaine from Hair

Sample Application

Apply the sample at a suitable flow rate (1-5ml/min).

Interference Elution

Rinse the column with water (2ml), hexane (2ml) and methanol (2ml) at a flow rate of 1-2ml/min. Dry the column prior to elution for 5 minutes.

Analyte Elution

Elute analytes with 78:20:2 v/v/v dichloromethane/isopropanol/ammonium hydroxide (2 x 1.5ml) at a flow rate of 1-2ml/min, including a 1 minute soak step after each aliquot.

Evaporation and Reconstitution

Evaporate under nitrogen.

Final Analysis

GC-MS with derivatisation.

Reagents

1. Water
2. Acetone
3. Ethanol
4. Methanol
5. 1.0M HCl
6. 0.1M phosphate buffer pH 6.0
7. Dichloromethane
8. Isopropanol
9. Ammonium hydroxide

UNITED KINGDOM

Kinesis Ltd

Tel: +44 (0)1480 212122
Fax: +44 (0)1480 212111
E-mail: sales@kinesis.co.uk
Web: kinesis.co.uk

GERMANY, SWITZERLAND & AUSTRIA

Kinesis GmbH (formerly Abimed)
Tel: +49 (0)2173 89 05-0
Fax: +49 (0)2173 89 05-77
Email: sales@kinesisgmbh.de
Web: kinesisgmbh.de

USA & CANADA

Kinesis Inc

Tel (toll free): (866) 934-6353
Tel: (518) 289-5817
Fax: (518) 289-5818
Email: sales@kinesis-usa.com
Web: kinesis-usa.com

AUSTRALIA & NEW ZEALAND

Kinesis Australia Pty Ltd

Tel: +61 (0)7 3829 3996
Fax: +61 (0)7 3829 3997
Email: sales@kinesis-australia.com.au
Web: kinesis-australia.com.au

