

8.7 OASIS

pipet4u accu pipet controller

Universal for all plastic or glass pipettes



High Comfort

The OASIS pipette controller is an engine powered pipetting aid designed for cordless work with glass or plastic pipettes in 0.1 – 100 ml volume range. Ergonomical shaped and carefully designed light weight handle, together with smooth push buttons and conveniently located slide switches ensure effortless pipetting even during extensive use.

Maximum Speed and Accuracy

It has a high performance micro pump with a powerful but quite engine with 3 speed modes. With soft grip-stress free control knobs it can aspirate-dispense 25 ml pipettes under 4 seconds. The pump speed can be fine-tuned by varying finger pressure to improve the sensitivity and for highest accuracy also in low volumes.

Safety System

It is equipped with replaceable membrane filter, which protect the device and the user from contamination by vapors of solutions being drawn in the pipette.

An incredible combination of ergonomic and value with all modern conveniences !

AHN OASIS pipette controller will be delivered complete with 2 x 1.2 V / 1000 mAh environmental friendly rechargeable NiMH batteries (cadmium free), with a battery charger, with 0.45 µm and 0.20 µm hydrophobic PTFE membrane filters and a wall bracket system for a safety handling. The intelligence charger prevents over charging and 2-3 hours to fully charge. The device has a low battery LED indicator. If the IED lamp lights up, it means that the batteries are discharged and they should be charged The OASIS has a fully autoclavable nosepiece and a continuous working limit of 10 hours !

Cat.no.	Description	Volume Range	Pack type	Sales unit
8-300-00-9	AHN OASIS Accu pipet controller	0,1-100 ml	carton	carton / 1 OASIS
8-300-03-0	AHN OASIS PTFE membrane filters	0.20 µm	carton	carton / 5 filters
8-300-04-0	AHN OASIS PTFE membrane filters	0.45 µm	carton	carton / 5 filters