


Hypersil GOLD Amino HPLC Columns


A high performance aminopropyl phase that gives excellent chromatographic properties in three modes: weak anion exchange, reversed phase and normal phase

Hypersil GOLD Amino columns have an aminopropyl ligand bonded to highly pure base deactivated silica.

- ▶ **Retains anions and organic acids in weak anion exchange**
- ▶ **Excellent for carbohydrate analysis in reversed phase**
- ▶ **Alternative selectivity to silica columns in normal phase chromatography**
- ▶ **Outstanding peak shape and sensitivity**
- ▶ **1.9µm particle size columns can be used to improve speed and efficiency**



Hypersil GOLD Amino HPLC Columns					
Length (mm)	1.0mm I.D.	2.1mm I.D.	3.0mm I.D.	4.0mm I.D.	4.6mm I.D.
Particle Size 1.9µm					
20mm	25702-021030	25702-022130	25702-023030	--	--
30mm	25702-031030	25702-032130	25702-033030	--	--
50mm	25702-051030	25702-052130	25702-053030	--	25702-054630
100mm	25702-101030	25702-102130	25702-103030	--	--
150mm	--	25702-152130	--	--	--
200mm	--	25702-202130	--	--	--
Particle Size 3µm					
30mm	25703-031030	25703-032130	25703-033030	25703-034030	25703-034630
50mm	25703-051030	25703-052130	25703-053030	25703-054030	25703-054630
100mm	25703-101030	25703-102130	25703-103030	25703-104030	25703-104630
150mm	25703-151030	25703-152130	25703-153030	25703-154030	25703-154630
Particle Size 5µm					
30mm	25705-031030	25705-032130	25705-033030	25705-034030	25705-034630
50mm	25705-051030	25705-052130	25705-053030	25705-054030	25705-054630
100mm	25705-101030	25705-102130	25705-103030	25705-104030	25705-104630
150mm	25705-151030	25705-152130	25705-153030	25705-154030	25705-154630
250mm	25705-251030	25705-252130	25705-253030	25705-254030	25705-254630

Hypersil GOLD Amino Drop-In Guard Cartridges						
Particle Size	Length	4.6mm/ 4.0mm I.D.	3.0mm I.D.	2.1mm I.D.	1.0mm I.D.	Quantity
3µm	10mm	25703-014001	25703-013001	25703-012101	25703-011001	4 Pack
5µm	10mm	25705-014001	25705-013001	25705-012101	25705-011001	4 Pack
	UNIGUARD Drop-in Guard Cartridge Holder	850-00	852-00	852-00	851-00	1 Each