


# Hypersil GOLD C<sub>4</sub> HPLC Columns


Lower hydrophobicity than C<sub>18</sub> or C<sub>8</sub> recommended for very hydrophobic analytes


Hypersil GOLD C<sub>4</sub> columns provide similar selectivity to C<sub>18</sub> and C<sub>8</sub> columns but with reduced retention.

- ▶ Lower hydrophobicity
- ▶ Faster separations
- ▶ Excellent peak shape
- ▶ High efficiency
- ▶ Outstanding sensitivity
- ▶ 1.9µm particle size columns can be used to improve speed and efficiency



Hypersil GOLD C <sub>4</sub> Analytical HPLC Columns					
Length	1.0mm I.D.	2.1mm I.D.	3.0mm I.D.	4.0mm I.D.	4.6mm I.D.
<b>Particle Size 1.9µm</b>					
20mm	25502-021030	25502-022130	25502-023030	--	--
30mm	25502-031030	25502-032130	25502-033030	--	--
50mm	25502-051030	25502-052130	25502-053030	--	25502-054630
100mm	25502-101030	25502-102130	25502-103030	--	--
150mm	--	25502-152130	--	--	--
200mm	--	25502-202130	--	--	--
<b>Particle Size 3µm</b>					
30mm	25503-031030	25503-032130	25503-033030	25503-034030	25503-034630
50mm	25503-051030	25503-052130	25503-053030	25503-054030	25503-054630
100mm	25503-101030	25503-102130	25503-103030	25503-104030	25503-104630
150mm	25503-151030	25503-152130	25503-153030	25503-154030	25503-154630
<b>Particle Size 5µm</b>					
30mm	25505-031030	25505-032130	25505-033030	25505-034030	25505-034630
50mm	25505-051030	25505-052130	25505-053030	25505-054030	25505-054630
100mm	25505-101030	25505-102130	25505-103030	25505-104030	25505-104630
150mm	25505-151030	25505-152130	25505-153030	25505-154030	25505-154630
250mm	25505-251030	25505-252130	25505-253030	25505-254030	25505-254630

Hypersil GOLD C <sub>4</sub> Drop-In Guard Cartridges						
Particle Size	Length	4.6mm/ 4.0mm I.D.	3.0mm I.D.	2.1mm I.D.	1.0mm I.D.	Quantity
3µm	10mm	25503-014001	25503-013001	25503-012101	25503-011001	4 Pack
5µm	10mm	25505-014001	25505-013001	25505-012101	25503-011001	4 Pack
	UNIGUARD Drop-in Guard Cartridge Holder	850-00	852-00	852-00	851-00	1 Each



**Go to: PAGE 90**

See our range of HyperSep SPE Products