VHG Labs brand Petroleum CRMs and reference materials



The VHG Labs brand oil standards have become the reference materials of choice in oil labs around the world.

We offer an extensive line of petroleum standards that are ideally suited for analysis by ICP, ICP-MS, XRF, UVF, and physical testing instruments.

Our products includes single and multi-element standards (both sulfonate-based and sulfur-free versions) as well as custom blends. Some of our most popular oil standards include:

- Sulfur Standards
 - Sulfur in a variety of matrices including: #2 diesel fuel, light mineral oil, heavy mineral oil, residual oil, isooctane, biodiesel and others. Suitable for use with many ASTM Methods
- Lube Oil Standards Suitable for ASTM D4927, D6481, D6443
- Chlorine in Oil Standards
 Suitable for ASTM D4929 or D5384

- Sulfur and Metals in Oil Standards Suitable for ASTM D5708
- Petroleum Physical Test Standards
 Viscosity, Flash Point, Pour Point, Cloud Point,
 Freezing Point, and Cold Filter Plug Point.
 Suitable for use with ASTM methods
- Biodiesel Standards
 Metals in biodiesel for ASTM D6751, EN 14107,
 EN 14108, EN 14109, EN 14214, EN 14538
- Proficiency Testing
 Looking to check the accuracy of your results or
 gain a competitive edge by demonstrating your
 competence? Ask about our proficiency testing
 programs for new and used oil.



For further information, or if you require substances or materials not currently listed please contact one of our local sales offices.

Brazil

Tel: +55 12 3302 5880 Email: bz@lgcstandards.com

Bulgaria

Tel: +359 (0)2 971 4955 Email: bg@lgcstandards.com

China

Tel: 800 810 4630 (toll-free) Email: infochina@lgcgroup.com

France

Tel: +33 (0)3 88 04 82 82 Email: fr@lgcstandards.com

Germany

Tel: +49 (0)281 9887 0 Email: de@lgcstandards.com

Hungary

Tel: +49 (0)281 9887 0 Email: de@lgcstandards.com

India

Tel: +91 (0)80 6701 2000 Email: in@lgcpromochem.com **Ireland**

Tel: +44 7879556983 Email: ire@lgcstandards.com

Italy

Tel: +39 02 2247 6412 Email: it@lgcstandards.com

Middle East

Tel: +49 (0)281 9887 0 Email: de@lgcstandards.com

Netherlands

Tel: +31 (0)643 775 422 Email: nl@lgcstandards.com

Nordic countries

Tel: +49 (0)281 9887 0 Email: de@lgcstandards.com

Poland

Tel: +48 22 751 31 40 Email: pl@lgcstandards.com

Romania

Tel: +40 364 116890

Email: ro@lgcstandards.com

Russia

Tel: +7 812 777 04 88 Email: ru@lgcgroup.com

South Africa

Tel: +27 (0)11 466 4321

Email: info@industrialanalytical.

co.za

Spain

Tel: +34 (0)93 308 4181 Email: es@lgcstandards.com

Turkey

Tel: +90 216 360 0870

Email: tur@lgcstandards.com

United Kingdom

Tel: +44 (0)20 8943 8480

Email: uksales@lgcstandards.com

USA + Canada

Tel: 1-855-LGC-USA1 (toll-free) Email: lgcusa@lgcstandards.com

www.lgcstandards.com • www.vhglabs.com





