

Detector Lamps & Instrument Spares

High Quality Deuterium UV Lamps & Pump Spares
for all Leading HPLC Instrumentation Brands



KINESIS LONG LIFE LAMPS FOR ALL LEADING HPLC UV DETECTORS INCLUDING:



Agilent

1100 1200 1260 1290

Kinesis Part Number	OEM Equivalent	Detector Models	Kinesis Guarantee	Matches OEM
LD-AGI-104LL	G1314-60100, G1314-60101	G1314A/B/D/E/F VWD	2000Hrs/1 Year	✓
LD-AGI-105LL	5181-1530 5182-1530 2140-0590	G1315A/B DAD G1365A/B DAD	2000Hrs/1 Year	✓
LD-AGI-107LL	2140-0820	G1315C/D DAD G1365C/D DAD	2000Hrs/1 Year	✓
LD-AGI-108LL	5190-0917	G4212A/B (8-Pin)	2000Hrs/1 Year	✓

Waters

996 2996 2487 Alliance

Kinesis Part Number	OEM Equivalent	Detector Models	Kinesis Guarantee	Matches OEM
LD-WAT-104LL	WAT052586	996 2996	2000Hrs/1 Year	✓
LD-WAT-105LL	WAS081142	2487 Alliance	2000Hrs/1 Year	✓

Shimadzu

SPD10A AVP AVVP M10AVP 20A 20AV 2010

Kinesis Part Number	OEM Equivalent	Detector Models	Kinesis Guarantee	Matches OEM
LD-SHI-101LL	228-34016-02	SPD10A AVP AVVP M10AVP 20A 20AV	2000Hrs/1 Year	✓
LD-SHI-103LL	228-37401-00	2010	2000Hrs/1 Year	✓

Dionex

Ultimate 3000 3000RS 3100 3400 PDA100 PDA3000
Gynkotek UVD 320 340S 160 170S 170U

Kinesis Part Number	OEM Equivalent	Detector Models	Kinesis Guarantee	Matches OEM
LD-DIO-102LL	939016T	PDA100 PDA3000	2000Hrs/1 Year	✓
LD-DIO-105LL	6074.1110	ULTIMATE	2000Hrs/1 Year	✓
LD-GYN-100LL	5053.1200	UVD 320 340S 160 170S 170U	2000Hrs/1 Year	✓

Merck Hitachi

LaChrom

Kinesis Part Number	OEM Equivalent	Detector Models	Kinesis Guarantee	Matches OEM
LD-MEH-100LL	892-2550	All L and U Series	2000Hrs/1 Year	✓

We have a comprehensive range of lamps not listed here all available on request for many other makes and models of UV HPLC detectors and Spectrophotometers.

Please contact Kinesis for full listing or visit kinesis-group.com

Guaranteed Quality Long Life HPLC Lamps

Deuterium 2000 Hour Lamps for HPLC UV Detection

Our Lamps

Kinesis have been manufacturing lamps since 2002 and we now have over 15 years experience in precision lamp alignment. Using processes verified by the world's leading lamp manufacturers, Kinesis are able to offer the highest quality lamp on the market today whilst offering customers huge price savings.

With sales last year of over 10,000 lamps this makes Kinesis the world's largest 3rd party lamp manufacturer with repeat business in every market sector across the globe.

Key Points

- 15 years lamp manufacturing experience
- World's largest 3rd party lamp producer
- Over 10,000 lamps sold every year
- Price saving benefit



Kinesis Quality Guarantee

By working closely with all the leading Deuterium lamp manufacturers, Kinesis are able to use in our alignment process the same bare lamp used by the original equipment manufacturer in the production of their branded detector lamps.

This enables Kinesis to offer, with complete confidence, the same warranties and guarantees you would receive when purchasing branded lamps from the original equipment manufacturers (OEM).

Key Points

- OEM equivalent lamps
- Same bare lamps as used to produce OEM lamp
- Guarantees and warranties match those offered by OEM

Certified Lamp Production

Kinesis lamps are manufactured using industry verified processes following stringently controlled ISO 9000 certified procedures. Every lamp we produce is traced by its serial number through the manufacturing process and tested, with results recorded, at every step – with a final QC intensity check in the original manufacturer's detector. Only at this point is the lamp certified and released for sale.

Key Points

- ISO 9000 certified manufacturing procedures
- Industry verified processes
- Lamps QC tested in original detector for final intensity readings
- Lamps certified for quality at every step



Kinesis HPLC Instrument Spares

High quality HPLC pump and autosampler spares for all leading instrument models including:

Agilent, ABI, Beckman, Dionex, Gilson, Jasco, Knauer, Merck Hitachi, Perkin Elmer, Shimadzu, Thermo, Varian and Waters



Pump Spares

- Piston seals
- Check valves
- Piston assemblies

Autosampler/Injection Spares

- Needles
- Syringes
- Needle seats
- Rotor seals

IDEX Accessories

- Tubing
- Fittings
- Filters

Maintenance Kits

- Kinesis quality preventative maintenance kits for instrumentation service including pump and autosampler spares
- Maintenance kits now available for the following instrumentation
 - Agilent
 - Shimadzu
 - Waters
- Other kits to follow

Price and Quality

- Highest quality parts
- Huge price savings



To ensure the highest product quality possible, Kinesis source their spare parts from the same manufacturers that are approved by many of the world's leading brands.

sales@kinesis-group.co.uk

UNITED KINGDOM

Kinesis Ltd

Tel: +44 (0)1480 212122

kinesis.co.uk

GERMANY, SWITZERLAND
& AUSTRIA

Kinesis GmbH

Tel: +49 (0)2173 89 05-0

kinesisgmbh.de

USA & CANADA

Kinesis Inc

Tel: (518) 289-5817

kinesis-usa.com

AUSTRALIA
& NEW ZEALAND

Kinesis Australia Pty Ltd

Tel: +61 (0)7 3829 3996

kinesis-australia.com.au

ITALY

ABreg S.r.l

Tel: +39 0131 389513

abreg.com



Kinesis recognise all trademarks and registered trademarks. KLTD-1803-0517

 **Kinesis**
Scientific Experts