

MarvelX™ Tubing Volumes

Material	ID in MM	Length in MM	Part Number	Description	Volume uL
PLS	0.025	70	UPFP-6025070	MarvelX PLS 25µmx70mm Kit	0.034
PLS	0.025	150	UPFP-6025150	MarvelX PLS 25µmx150mm Kit	0.074
PLS	0.025	250	UPFP-6025250	MarvelX PLS 25µmx250mm Kit	0.123
PLS	0.025	350	UPFP-6025350	MarvelX PLS 25µmx350mm Kit	0.172
PLS	0.025	500	UPFP-6025500	MarvelX PLS 25µmx500mm Kit	0.246
PLS	0.025	600	UPFP-6025600	MarvelX PLS 25µmx600mm Kit	0.295
PLS	0.050	70	UPFP-6050070	MarvelX PLS 50µmx70mm Kit	0.137
PLS	0.050	150	UPFP-6050150	MarvelX PLS 50µmx150mm Kit	0.295
PLS	0.050	250	UPFP-6050250	MarvelX PLS 50µmx250mm Kit	0.491
PLS	0.050	350	UPFP-6050350	MarvelX PLS 50µmx350mm Kit	0.687
PLS	0.050	500	UPFP-6050500	MarvelX PLS 50µmx500mm Kit	0.982
PLS	0.050	600	UPFP-6050600	MarvelX PLS 50µmx600mm Kit	1.178
PLS	0.075	70	UPFP-6075070	MarvelX PLS 75µmx70mm Kit	0.309
PLS	0.075	150	UPFP-6075150	MarvelX PLS 75µmx150mm Kit	0.663
PLS	0.075	250	UPFP-6075250	MarvelX PLS 75µmx250mm Kit	1.105
PLS	0.075	350	UPFP-6075350	MarvelX PLS 75µmx350mm Kit	1.547
PLS	0.075	500	UPFP-6075500	MarvelX PLS 75µmx500mm Kit	2.210
PLS	0.075	600	UPFP-6075600	MarvelX PLS 75µmx600mm Kit	2.652
PLS	0.100	70	UPFP-6100070	MarvelX PLS 100µmx70mm Kit	0.550
PLS	0.100	150	UPFP-6100150	MarvelX PLS 100µmx150mm Kit	1.178
PLS	0.100	250	UPFP-6100250	MarvelX PLS 100µmx250mm Kit	1.964
PLS	0.100	350	UPFP-6100350	MarvelX PLS 100µmx350mm Kit	2.750
PLS	0.100	500	UPFP-6100500	MarvelX PLS 100µmx500mm Kit	3.928
PLS	0.100	600	UPFP-6100600	MarvelX PLS 100µmx600mm Kit	4.714
PLS	0.150	70	UPFP-6150070	MarvelX PLS 150µmx70mm Kit	1.237
PLS	0.150	150	UPFP-6150150	MarvelX PLS 150µmx150mm Kit	2.652
PLS	0.150	250	UPFP-6150250	MarvelX PLS 150µmx250mm Kit	4.419
PLS	0.150	350	UPFP-6150350	MarvelX PLS 150µmx350mm Kit	6.187
PLS	0.150	500	UPFP-6150500	MarvelX PLS 150µmx500mm Kit	8.839
PLS	0.150	600	UPFP-6150600	MarvelX PLS 150µmx600mm Kit	10.606
PLS	0.300	70	UPFP-6300070	MarvelX PLS 300µmx70mm Kit	4.950
PLS	0.300	150	UPFP-6300150	MarvelX PLS 300µmx150mm Kit	10.606
PLS	0.300	250	UPFP-6300250	MarvelX PLS 300µmx250mm Kit	17.677
PLS	0.300	350	UPFP-6300350	MarvelX PLS 300µmx350mm Kit	24.748
PLS	0.300	500	UPFP-6300500	MarvelX PLS 300µmx500mm Kit	35.354
PLS	0.300	600	UPFP-6300600	MarvelX PLS 300µmx600mm Kit	42.425
SS	0.100	70	UPFS-6100070	MarvelX SS 100µmx70mm Kit	0.550
SS	0.100	150	UPFS-6100150	MarvelX SS 100µmx150mm Kit	1.178
SS	0.100	250	UPFS-6100250	MarvelX SS 100µmx250mm Kit	1.964
SS	0.100	350	UPFS-6100350	MarvelX SS 100µmx350mm Kit	2.750
SS	0.100	500	UPFS-6100500	MarvelX SS 100µmx500mm Kit	3.928
SS	0.100	600	UPFS-6100600	MarvelX SS 100µmx600mm Kit	4.714
SS	0.125	70	UPFS-6125070	MarvelX SS 125µmx70mm Kit	0.859
SS	0.125	150	UPFS-6125150	MarvelX SS 125µmx150mm Kit	1.841
SS	0.125	250	UPFS-6125250	MarvelX SS 125µmx250mm Kit	3.069
SS	0.125	350	UPFS-6125350	MarvelX SS 125µmx350mm Kit	4.297
SS	0.125	500	UPFS-6125500	MarvelX SS 125µmx500mm Kit	6.138
SS	0.125	600	UPFS-6125600	MarvelX SS 125µmx600mm Kit	7.365
SS	0.254	70	UPFS-6254070	MarvelX SS 254µmx70mm Kit	3.548
SS	0.254	150	UPFS-6254150	MarvelX SS 254µmx150mm Kit	7.603
SS	0.254	250	UPFS-6254250	MarvelX SS 254µmx250mm Kit	12.672
SS	0.254	350	UPFS-6254350	MarvelX SS 254µmx350mm Kit	17.740
SS	0.254	500	UPFS-6254500	MarvelX SS 254µmx500mm Kit	25.343
SS	0.254	600	UPFS-6254600	MarvelX SS 254µmx600mm Kit	30.412
PLS	0.025	70	UPFP-6025070T	MarvelX PLS 25µmx70mm Tube	0.034
PLS	0.025	150	UPFP-6025150T	MarvelX PLS 25µmx150mm Tube	0.074
PLS	0.025	250	UPFP-6025250T	MarvelX PLS 25µmx250mm Tube	0.123
PLS	0.025	350	UPFP-6025350T	MarvelX PLS 25µmx350mm Tube	0.172
PLS	0.025	500	UPFP-6025500T	MarvelX PLS 25µmx500mm Tube	0.246
PLS	0.025	600	UPFP-6025600T	MarvelX PLS 25µmx600mm Tube	0.295
PLS	0.050	70	UPFP-6050070T	MarvelX PLS 50µmx70mm Tube	0.137
PLS	0.050	150	UPFP-6050150T	MarvelX PLS 50µmx150mm Tube	0.295
PLS	0.050	250	UPFP-6050250T	MarvelX PLS 50µmx250mm Tube	0.491
PLS	0.050	350	UPFP-6050350T	MarvelX PLS 50µmx350mm Tube	0.687
PLS	0.050	500	UPFP-6050500T	MarvelX PLS 50µmx500mm Tube	0.982
PLS	0.050	600	UPFP-6050600T	MarvelX PLS 50µmx600mm Tube	1.178
PLS	0.075	70	UPFP-6075070T	MarvelX PLS 75µmx70mm Tube	0.309
PLS	0.075	150	UPFP-6075150T	MarvelX PLS 75µmx150mm Tube	0.663
PLS	0.075	250	UPFP-6075250T	MarvelX PLS 75µmx250mm Tube	1.105
PLS	0.075	350	UPFP-6075350T	MarvelX PLS 75µmx350mm Tube	1.547
PLS	0.075	500	UPFP-6075500T	MarvelX PLS 75µmx500mm Tube	2.210
PLS	0.075	600	UPFP-6075600T	MarvelX PLS 75µmx600mm Tube	2.652
PLS	0.100	70	UPFP-6100070T	MarvelX PLS 100µmx70mm Tube	0.550
PLS	0.100	150	UPFP-6100150T	MarvelX PLS 100µmx150mm Tube	1.178
PLS	0.100	250	UPFP-6100250T	MarvelX PLS 100µmx250mm Tube	1.964
PLS	0.100	350	UPFP-6100350T	MarvelX PLS 100µmx350mm Tube	2.750
PLS	0.100	500	UPFP-6100500T	MarvelX PLS 100µmx500mm Tube	3.928
PLS	0.100	600	UPFP-6100600T	MarvelX PLS 100µmx600mm Tube	4.714
PLS	0.150	70	UPFP-6150070T	MarvelX PLS 150µmx70mm Tube	1.237
PLS	0.150	150	UPFP-6150150T	MarvelX PLS 150µmx150mm Tube	2.652

Material	ID in MM	Length in MM	Part Number	Description	Volume uL
PLS	0.150	250	UPFP-6150250T	MarvelX PLS 150µmx250mm Tube	4.419
PLS	0.150	350	UPFP-6150350T	MarvelX PLS 150µmx350mm Tube	6.187
PLS	0.150	500	UPFP-6150500T	MarvelX PLS 150µmx500mm Tube	8.839
PLS	0.150	600	UPFP-6150600T	MarvelX PLS 150µmx600mm Tube	10.606
PLS	0.300	70	UPFP-6300070T	MarvelX PLS 300µmx70mm Tube	4.950
PLS	0.300	150	UPFP-6300150T	MarvelX PLS 300µmx150mm Tube	10.606
PLS	0.300	250	UPFP-6300250T	MarvelX PLS 300µmx250mm Tube	17.677
PLS	0.300	350	UPFP-6300350T	MarvelX PLS 300µmx350mm Tube	24.748
PLS	0.300	500	UPFP-6300500T	MarvelX PLS 300µmx500mm Tube	35.354
PLS	0.300	600	UPFP-6300600T	MarvelX PLS 300µmx600mm Tube	42.425
SS	0.100	70	UPFS-6100070T	MarvelX SS 100µmx70mm Tube	0.550
SS	0.100	150	UPFS-6100150T	MarvelX SS 100µmx150mm Tube	1.178
SS	0.100	250	UPFS-6100250T	MarvelX SS 100µmx250mm Tube	1.964
SS	0.100	350	UPFS-6100350T	MarvelX SS 100µmx350mm Tube	2.750
SS	0.100	500	UPFS-6100500T	MarvelX SS 100µmx500mm Tube	3.928
SS	0.100	600	UPFS-6100600T	MarvelX SS 100µmx600mm Tube	4.714
SS	0.125	70	UPFS-6125070T	MarvelX SS 125µmx70mm Tube	0.859
SS	0.125	150	UPFS-6125150T	MarvelX SS 125µmx150mm Tube	1.841
SS	0.125	250	UPFS-6125250T	MarvelX SS 125µmx250mm Tube	3.069
SS	0.125	350	UPFS-6125350T	MarvelX SS 125µmx350mm Tube	4.297
SS	0.125	500	UPFS-6125500T	MarvelX SS 125µmx500mm Tube	6.138
SS	0.125	600	UPFS-6125600T	MarvelX SS 125µmx600mm Tube	7.365
SS	0.254	70	UPFS-6254070T	MarvelX SS 254µmx70mm Tube	3.548
SS	0.254	150	UPFS-6254150T	MarvelX SS 254µmx150mm Tube	7.603
SS	0.254	250	UPFS-6254250T	MarvelX SS 254µmx250mm Tube	12.672
SS	0.254	350	UPFS-6254350T	MarvelX SS 254µmx350mm Tube	17.740
SS	0.254	500	UPFS-6254500T	MarvelX SS 254µmx500mm Tube	25.343
SS	0.254	600	UPFS-6254600T	MarvelX SS 254µmx600mm Tube	30.412
SS			UPN-61032	Replacement Fittings 3/Pack	