



FLAME Filler

Ergonomic fast pipetting solution with accurate speed control



www.gvs.com



Features

- ◆ Intuitive and convenient with 6 speed levels
- ◆ LCD screen displays battery status and pipetting speed
- ◆ Lightweight pipette filler compatible with glass or plastic serological pipettes (0.1 -100 mL)
- ◆ Quiet and powerful motor and pump fills 25 mL pipette in less than 5 seconds
- ◆ Ergonomic design for fatigue-free and effective streamlined pipetting
- ◆ Long-lasting lithium-ion battery

Specifications

| | |
|----------------------|--|
| Aspiration Speed | 25mL<5s (6 Level) |
| Dispensing Speed | Motor/Gravity |
| Battery | Replaceable lithium battery |
| Battery Service Life | Over 12 hours of intermittent use |
| Charging Time | 4-5 Hours |
| Pipette Types | Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes |
| Filter | 0.45µm Hydrophobic |
| Weight | 208g |