

Universal connection for (U)HPLC Columns

Connect with confidence using the EasyConn kit. No more under- or overtightening, no more uncertainty. Just tighten until you hear the click - and know instantly that the connection is perfectly assured.

Key Benefits of EasyConn

- **Guaranteed Perfect Tightening** - Eliminates the risk of under or overtightening. The audible click confirms a precise, secure connection every time.
- **Reliable Performance up to UHPLC Pressures** - Designed to withstand pressures up to 1300 bar, ensuring dependable operation even in demanding UHPLC applications.
- **Fast, Simple Installation** - Quick, tool free assembly with zero dead volume - just tighten until it clicks and you're ready to run.
- **Universal UHPLC/HPLC Compatibility** - Compatible with all common UHPLC and HPLC column ports (M4, M2.5, 1/16"), offering true system flexibility.
- **Robust and Durable Design** - Built to minimize the risk of damage during frequent handling and reconnections, extending component lifetime.
- **Compact Form Factor** - Small and space saving design, ideal for tight instrument layouts and crowded column compartments.
- **Excellent Chemical Resistance** - Available in inert materials, including bio compatible PEEK options for sensitive analytical applications.
- **Leak Free and Maintenance Free Operation** - No leaks, no complex adjustments - just a consistently secure connection you can trust.




EasyConn

Universal connection for (U)HPLC Columns

No leaks. No tedious adjustments.
Just fast, reliable, and safe connections.
EasyConn provides a simple, robust solution
for capillary-to-column connections, ensuring consistent
performance even at the highest pressures.

Product number		Availability
ECON-7254150	EasyConn SS Kit, 254 μ m x 150mm	in stock
ECON-7125600	EasyConn SS Kit, 125 μ m x 600mm	in stock
ECON-7125500	EasyConn SS Kit, 125 μ m x 500mm	in stock
ECON-7125350	EasyConn SS Kit, 125 μ m x 350mm	in stock
ECON-7125250	EasyConn SS Kit, 125 μ m x 250mm	in stock
ECON-7125150	EasyConn SS Kit, 125 μ m x 150mm	in stock
ECON-7100350	EasyConn SS Kit, 100 μ m x 350mm	in stock
ECON-7100250	EasyConn SS Kit, 100 μ m x 250mm	in stock
ECON-7100150	EasyConn SS Kit, 100 μ m x 150mm	in stock



Ask for your quote by email: j.suli@chromservis.eu or inquiries@chromservis.eu
For more details please log on our web: www.chromservis.eu